

# Total Protein

## Interpretive Summary

---

**Description:** Total protein in the blood is primarily composed of albumin and globulins (see individual interpretive summaries) with a minor contribution from coagulation proteins and fibrinogen.

---

### Decreased Total Protein

#### Common Causes

- Decreased production:
  - Intestinal malabsorption
  - Malnutrition
  - Severe gastrointestinal (GI) parasitism (multifactorial)
  - Liver insufficiency/failure
- Increased loss
  - Hemorrhage (especially external)
  - Protein-losing nephropathy (PLN)
  - Protein-losing enteropathy (PLE)/lymphangiectasia
- Other causes
  - Exocrine pancreatic insufficiency (EPI) leading to maldigestion

#### Uncommon Causes

- Decreased production
  - Addison's disease (decreased albumin production)
  - Immunodeficiency (decreased globulins)
    - Congenital
    - Inherited
    - Acquired
- Increased loss
  - Severe exudative skin disease
  - Vasculitis
  - Peritonitis
  - Pleuritis
- Other causes
  - Failure of passive transfer of colostrum (neonates)
  - Hemodilution
    - Intravenous fluids
    - Edema disorders
      - Congestive heart failure
      - Nephrotic syndrome
      - Hydrothorax
      - Ascites
    - Concurrent hypovolemia and increased total body water
      - Fluid accumulation in third space
      - Fluid accumulation in GI tract
    - Syndrome of inappropriate ADH secretion (SIADH)

#### Related Findings

- Decreased production
  - Intestinal malabsorption
    - Decreased albumin and globulins

- Albumin:globulin (A:G) ratio should be normal
    - Increased or decreased serum folate, decreased cobalamin
  - GI parasitism
    - Positive fecal ova and parasites
  - Liver disease
    - Increased liver enzymes (ALT, AST, GGT, ALP), bile acids
    - Decreased albumin, BUN, cholesterol, glucose
    - Primarily due to decreased albumin production; A:G ratio decreased
- Increased loss
  - Hemorrhage
    - Regenerative anemia on CBC
    - Decreased albumin and globulins
    - A:G ratio should be normal
  - Protein-losing nephropathy
    - Increased urine protein:creatinine ratio
    - Decreased urine specific gravity
    - Increased BUN, creatinine, phosphorus with secondary tubular damage
    - Primarily loss of albumin. A:G ratio decreased
    - Globulins may be increased if chronic antigenic stimulation due to underlying cause
  - Protein losing enteropathy
    - Both albumin and globulins are decreased
    - A:G ratio should be normal
    - Increased or decreased serum folate, decreased cobalamin (B12)
    - Fecal alpha one protease inhibitor increased
    - Decreased cholesterol with lymphangiectasia and also with other severe intestinal disease
- EPI
  - Decreased trypsin-like immunoreactivity (TLI)

## Increased Total Protein

### Common Causes

- Dehydration
- Chronic inflammation
- Infection
  - Viral (feline infectious peritonitis [FIP])
  - Rickettsial (especially ehrlichiosis)
  - Chronic bacterial
- Neoplasia
  - Plasma cell myeloma
  - Extra medullary plasmacytoma
  - Some lymphomas

### Uncommon Causes

- Infection
  - Viral (feline immunodeficiency virus [FIV])
  - Protozoal
  - Fungal
  - Parasites (e.g. heartworms)
- Immune-mediated disease
- Artifact (all cause false increases)
  - Lipemia
  - Hemolysis
  - Icterus

## Related Findings

- Dehydration
  - Increased hematocrit
  - Pre-renal azotemia
  - Well concentrated urine
  - Albumin and globulin both increased; A:G ratio normal
- Chronic inflammation
  - Increased neutrophils and monocytes
- Infection
  - Polyclonal gammopathy on serum protein electrophoresis (SPE); occasionally monoclonal in ehrlichiosis
  - Positive viral antibody or antigen titers/PCR
- Neoplasia
  - Monoclonal gammopathy on SPE
  - Bence Jones protein positive in urine (myeloma, some extra medullary plasma cell tumors)
  - Monoclonal gammopathy on urine protein electrophoresis (myeloma, some extra medullary plasma cell tumors)
  - CBC may show cytopenias in other cell lines
- Artifact
  - Gross examination of sample shows
    - Hemolysis
    - Icterus
    - Lipemia

---

## Additional Information

---

### Physiology

- Total protein in serum is made up of albumin and globulins only. In plasma there is also some contribution from fibrinogen and coagulation proteins
- Total protein in the body is manufactured mostly in the liver and by the structures that comprise the immune system (which make gamma-globulins)

### References

- Kaneko JJ, Harvey JW, Bruss ML. *Clinical Biochemistry of Domestic Animals*, 6th ed. San Diego, CA: Academic Press; 2008.
- Latimer KS, Mahaffey EA, Prasse KW, eds. *Duncan and Prasse's Veterinary Laboratory Medicine: Clinical Pathology*, 4th ed. Ames, IA: Blackwell; 2003.
- Stockham SL, Scott MA. *Fundamentals of Veterinary Clinical Pathology*, 2nd ed. Ames, IA: Blackwell; 2008.
- Willard MD, Tvedten H, eds. *Small Animal Clinical Diagnosis by Laboratory Methods*, 4th ed. St. Louis, MO: Saunders; 2004.

Last updated 11/1/2013