Diagnostic Update

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## Literaturnachweise

1. Miró G, Wright I, Michael H, et al. Seropositivity of main vector-borne pathogens in dogs across Europe. *Parasit Vectors*. 2022;15(1):189. doi:10.1186/s13071-022-05316-5
2. Stillman BA, Monn M, Liu J, et al. Performance of a commercially available in-clinic ELISA for detection of antibodies against *Anaplasma phagocytophilum*, *Anaplasma platys*, *Borrelia burgdorferi*, *Ehrlichia canis*, and *Ehrlichia ewingii* and *Dirofilaria immitis* antigen in dogs. *JAVMA*. 2014;245(1):80–86. doi:10.2460/javma.245.1.80
3. Chandrashekar R, Mainville CA, Beall MJ, et al. Performance of a commercially available in-clinic ELISA for the detection of antibodies against *Anaplasma phagocytophilum*, *Ehrlichia canis*, and *Borrelia burgdorferi* and *Dirofilaria immitis* antigen in dogs. *Am J Vet Res*. 2010;71(12):1443–1450. doi:10.2460/ajvr.71.12.1443
4. Beall M, Mainville C, Arguello-Marin A, et al. An improved point-of-care ELISA for the diagnosis of anaplasmosis and ehrlichiosis during the acute phase of tick-borne infections in dogs. *Topics in Companion Anim Med*. 2022;51:100735. doi:10.1016/j.tcam.2022.100735
5. Chandrashekar R, Beall MJ, Thatcher B, Saucier JM, Tyrrell P, Lappin MR. Serologic responses to peptides of *Anaplasma phagocytophilum* and *Borrelia burgdorferi* in dogs infested with wild-caught *Ixodes scapularis*. *Vet J*. 2017;226:6–11. doi:10.1016/j.tvjl.2017.06.005
6. Kova evi Filipovi MM, Beleti AD, Ili Božovi AV, et al. Molecular and serological prevalence of *Anaplasma phagocytophilum*, *A. platys*, *Ehrlichia canis*, *E. chaffeenses*, *E. ewingii*, *Borrelia burgdorferi*, *Babesia canis*, *B. gibsoni* and *B. vogeli* among clinically healthy outdoor dogs in Serbia. *Vet Parasitol Reg Stud Reports*. 2018;14:117–122. doi:10.1016/j.vprsr.2018.10.001
7. Daten hinterlegt bei IDEXX Laboratories, Inc. Westbrook, Maine USA.
8. Pantchev N, Pluta S, Huisinga E, et al. Tick-borne diseases (borreliosis, anaplasmosis, babesiosis) in German and Austrian dogs: status quo and review of distribution, transmission, clinical findings, diagnostics and prophylaxis. *Parasitol Res*. 2015;114 Suppl 1:S19–S54.
9. Gaunt S, Beall M, Stillman B, et al. Experimental infection and co-infection of dogs with *Anaplasma platys* and *Ehrlichia canis*: hematologic, serologic and molecular findings. *Parasit Vectors*. 2010;3(1):33. doi:10.1186/1756-3305-3-33
10. Scorpio DG, Dumler JS, Barat NC, et al. Comparative strain analysis of *Anaplasma phagocytophilum* infection and clinical outcomes in a canine model of granulocytic anaplasmosis. *Vector Borne Zoonotic Dis*. 2011;11(3):223–229. doi:10.1089/vbz.2009.0262
11. Richardson SS, Mainville CA, Arguello-Marin A, et al. A second-generation, point-of-care immunoassay provided improved detection of *Anaplasma* and *Ehrlichia* antibodies in PCR-positive dogs naturally infected with *Anaplasma* or *Ehrlichia* species. *J Vet Diagn Invest*. 2023;35(4):367–374. doi:10.1177/10406387231172723
12. Canine leishmaniosis: clinical staging, treatment and prognosis [fact sheet]. LeishVet website. Zugriff am 17.10.2023. [www.leishvet.org/fact-sheet/clinical-staging](http://www.leishvet.org/fact-sheet/clinical-staging)

### IDEXX GmbH

Humboldtstr. 2  
70806 Kornwestheim  
idexx.de

### IDEXX Diavet AG

Schlyffstr. 10  
8806 Bäch  
idexx.ch

### IDEXX Vet Med Labor GmbH

Börsegasse 12/1  
1010 Wien  
idexx.at